



Complex assessment of oxygen metabolism kinetics

The assessment of oxygen metabolism kinetics can be estimated by a PO_2 polarography technique (Chizhov 1983, 1987) on a PO_2 monitor.



A regional limb ischemia test is taken to measure oxygen consumption.

Polarographic curve in the assessment of oxygen metabolism kinetics



Measurement results

The program automatically generates the following indices

index
TcPo₂ , mm Hg
КСПК (OARC) , sec ⁻¹ (O ₂ absorption rate constant)
КСБК (ORRC) , sec ⁻¹ (O ₂ recovery rate constant)
ККР (COR) = (КСБК ORRC / КСПК OARC) (coefficient of oxygen reserve)
ККК (OCC) , mm Hg (O ₂ critical concentration)
t_{resp} , seconds (O ₂ response time)
ВАП, (АРТ) , sec (aerobic processes time)
ВАНП (AnPT) , sec (anaerobic processes time)
КААнГ (CAnGA) = (ВАНП AnPT / ВАП АРТ), sec (coefficient of anaerobic glycolysis activity)

The program also estimates the meaning of the indices presented in the table above and generates a conclusion characterising:

- oxygen microcirculation in tissues;
- tissue respiration processes;
- the functional backup of oxygen dependent enzyme systems;
- tissue oxygen consumption and delivery